

Glasschem reactor setup

Overhead stirrer mounted directly on flange lid. This gives the advantage of smooth stirring. Either stainless steel bearings or a magnetic coupled bearing for vacuum work.

Hot water/oil reservoir with internal pump. This system can pump hot oil up to 250 deg C through the jacketed glass vessels jacket. Precise temperature control via microprocessor controlled PID.



Control instruments according to customers specifications, usually temperature control, stirring speed, Ph, ORP or conductivity.

Jacketed glass vessel from 250ml up to 5000ml. Stainless steel or glass flanged lid to accommodate either temperature, Ph, ORP or any other probes. Sampling ports as per request

